|            |            |            |
|------------|------------|------------|
| <b>REF</b> |            |            |
| 05-1000/01 | 05-1001/03 | 05-1002/03 |
| 05-1003/03 | 05-1030    | 05-2001/03 |
| 05-2002/03 | 05-2003/03 | 05-2030    |
| <b>CE</b>  |            |            |

DEUTSCH

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#### EL

#### EN

This document does not contain all of the information necessary for application of the instruments. For proper handling, refer to other device-related instructions, e.g., the associated surgical technique, other instructions for use, such as the instructions for use H50, device-specific handling recommendations, and cleaning instructions as well as device labels.

Refer to the identification tag on the instrument and/or the packaging label for definitive identification information for the device, e.g., system compatibility, article number, materials, and shelf life. You should also take advantage of the training courses and printed materials provided for your information. To learn more, please contact the Waldemar Link GmbH & Co. KG sales office or your field representative.

#### Intended Use

The sterilization containers are used as re-usable packaging for sterilization with moist heat with medical devices. The sterilization container is made out of metal and complies with EN ISO 11607-1 and EN 868-8. To ensure contact between the items being sterilized and the sterilization medium, the sterilization container is provided with a single use filter made of paper. The sterilization container is designed for sterilization with moist heat as defined in DIN EN ISO 17665 and is suitable for use in sterilizers as defined in DIN EN 285 or DIN EN 58946. Please follow the instructions of your sterilizer manufacturer.

The quality of the sterility of the sterilized items depends in part on the storage conditions that prevail at the operator's site and must be specified with the hygiene officer of the operator. We recommend taking into account the environmental conditions defined in DIN 58953-9 for the storage of the sterilization containers.

- Warnings**
  - The sterilization containers are designed for sterilization with moist heat (steam sterilization) using the fractionated pre-vacuum method. The temperature must not exceed 137 °C because otherwise handles, insulation, or other non-metallic components may be damaged.
  - To ensure that the lid can easily close, the sterilization container should only be loaded up to max. 20 mm below the edge of the container.
  - For reasons of sterilization safety, the total weight of each sterilization container must not exceed 10 kg.
  - If prior contamination is suspected, the recommendations of the Robert Koch Institute (www.rki.de) must be followed. Further information on the proper reprocessing of medical devices can also be obtained from the "Arbeitskreis Instrumenten-Aufbereitung" [Working Committee Instrument Reprocessing] at www.a-k-i.org, specifically, the red brochure "Instrument reprocessing done correctly" and the gray brochure "Test series and statements".
  - For information on the correct use to Waldemar LINK GmbH & Co. KG for servicing purposes must be reprocessed beforehand by the sender since there is no hazard for third parties resulting from the devices.

**First Use**
The sterilization containers must be thoroughly cleaned, disinfected, and sterilized before the first use.

- Application**
  - Only sterilization containers without visible damage may be used. This applies particularly to the sealing surface of the container seal on the lid and to the container edge. Likewise, bent lids or filter mats must be immediately replaced because otherwise sterile storage of the sterilized items is not assured.
  - Before each use, it must be ensured that there is no solid or liquid contamination in or on the product. The functionality of the filter must be checked (visual inspection). The filter must not be damaged and should sit flush in the lid.
  - All clasps on the sterilization container must be checked for a firm seat.
  - The tension locks on the sides of the sterilization container must close with noticeable tension when the lid lies fully on top so that the microbial barrier of the container is assured.
  - We recommend additionally wrapping the loaded tray inserts with suitable cloths or liners for the sterilization.

- Medical devices that are sent to Waldemar LINK GmbH & Co. KG for servicing purposes must be reprocessed beforehand by the sender since there is no hazard for third parties resulting from the devices.
- For reasons of sterilization safety, the total weight of each sterilization container must not exceed 10 kg.
- If prior contamination is suspected, the recommendations of the Robert Koch Institute (www.rki.de) must be followed. Further information on the proper reprocessing of medical devices can also be obtained from the "Arbeitskreis Instrumenten-Aufbereitung" [Working Committee Instrument Reprocessing] at www.a-k-i.org, specifically, the red brochure "Instrument reprocessing done correctly" and the gray brochure "Test series and statements".
- For information on the correct use to Waldemar LINK GmbH & Co. KG for servicing purposes must be reprocessed beforehand by the sender since there is no hazard for third parties resulting from the devices.

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