	DE	Gebrauchsanweisung MobileLink Hüftpfannensystem
	EN	Instructions for use MobileLink Acetabular Cup System
	BG	Инструкции за употреба Системата с ацетабуларна чашка MobileLink
	CS	Návod k použití Systém acetabulární jamky MobileLink
	DA	Bruksanvisning MobileLink Acetabular Cup System
	EL	Οδηγίες χρήσης σύστημα κοτυλιαίου κυπελλίου MobileLink
	ES	Instrucciones de uso Sistema de cotilo acetabular MobileLink
	ET	Kasutusjuhend Puusanapasüsteem MobileLink
	FI	Käyttöohje MobileLink-lonkkamaljakuppijärjestelmä
	FR	Mode d'emploi Système de cupule acétabulaire MobileLink
	HR	Upute za korištenje Sustav MobileLink Acetabular Cup System
	HU	Használati útmutató MobileLink Acetabular Cup System
	ID	Instruksi Penggunaan Sistem Mangkuk MobileLink Acetabular
	IT	Istruzioni per l'uso Sistema di coppa acetabolare MobileLink
	LT	Naudojimo instrukcijos MobileLink gūžduobės taurelės sistema
	LV	Lietošanas instrukcija MobileLink acetabulārā komponenta sistēma
	NL	Gebruiksaanwijzing MobileLink acetabulair cupsysteem
	NO	Bruksanvisning MobileLink-leddskålsystemet
	PL	Instrukcja użycia System czaszy panewki stawowej MobileLink
	PT	Instruções de utilização Sistema de taça acetabular MobileLink
	RO	Instructiuni de folosire Sistemul de cupă acetabulară MobileLink
	RU	Инструкция по применению Имплантаты Система ацетабулярных чашек MobileLink
	SK	Inštrukcie na používanie Systém acetabulárnej jamky MobileLink
	SL	Navodila za uporabo Sistem acetabularne skodelice MobileLink
	SV	Bruksanvisning Implantat MobileLink höftledsskålsystem
	TR	Kullanım Talimatı İmplantlar MobileLink Asetabular Kap Sistemi
	UK	Інструкція для застосування Система вертлюжної чаші MobileLink
	ZH	说明书 骨科植入物 MobileLink 髌臼杯系统

# LINK®

CE<sub>0482</sub> / CE

Not for US distribution.

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## 1. General

Please read this document carefully before using the system and keep it for future consultation! Failure to comply with these instructions for use may cause damage to patient's health, implants and instruments.

This document does not contain all of the information necessary for the selection and application of the system. For safe and proper handling, refer to further product-related instructions, such as the surgical technique associated with the system, the packaging and cleaning instructions and the device labels on the packaging. For instruments that require calibration, please refer to the Additional Product Information, Torque Wrench (IFU-EU-134-140-001).

## 2. User Group and Environment

The intended users of this system are experienced and trained

- surgeons,
- theatre nurses,
- CSSD staff (CSSD = Central Sterile Supply Department).

## 3. Patient Group

The patient group for this medical device consists of adult, anaesthetised patients of any ethnic origin and of any gender, in whom one or more of the described indications are present and who are not considered unsuitable on the grounds of the listed contraindications.

## 4. System Description

The MobileLink Acetabular Cup System is part of a mechanical reconstruction of the hip joint.

The MobileLink Acetabular Cup System is an uncemented press-fit acetabular cup system, that consists of different components.

The base are press-fit shell components, that are fixed in the acetabulum. Into these shells different liners are to be introduced to give the intersection to the femoral head components. In between the shell and liner metal adapters – Shell/Insert Adapters - can be positioned to change the inclination of the cup entry plane and/or offset of the centre of rotation. With a special Dual Mobility Insert the MobileLink Acetabular Cup System can be transformed into a modular Dual Mobility System.

The instruments are designed for the implantation of the MobileLink Acetabular Cup System. The instruments are bundles on different instrument trays:

- Basic Instruments
- Revision Instruments
- Instruments for Size 74 – 80 mm
- Instruments for Dual Mobility Inserts (Option 1, Option 2)

In addition, the surgeon is supplied with the complementary acetabular reamers for the preparation of the acetabulum. Furthermore, a Trial Head instrument set is supplied on request. The provided instrument sets for the implantation of all different components of the MobileLink Acetabular Cup System enable the surgeon to perform the preparation of the acetabulum, a trial reposition as well as the final implantation of the system.

Next to the instrument sets, a selection of additional, optional instruments is available upon request.

## 5. Intended Use

### Implants

The non-active, surgically-invasive implantable MobileLink Acetabular Cup System manufactured by Waldemar Link GmbH & Co. KG is intended for long-term replacement of the acetabular side of a diseased and / or defective hip joint in the human body. The MobileLink Acetabular Cup System forms a total replacement of the hip joint when combined with the prosthesis head and stem. The MobileLink Acetabular Cup System can be used with full-grown, anaesthetised patients of any ethnic origin and sex. The MobileLink Acetabular Cup System is implanted without cement.

The implants may only be used and operated in an aseptic medical environment by persons who have the required training, knowledge and experience in the orthopedic and surgical field. The implants are supplied in sterile condition individually packed as single-use products.

### Instruments

The purpose of the instrumentation is to allow the user to use the associated implant system as part of the procedures described in the associated surgical technique. Any other use of the instrumentation is not permitted.

The instruments consist of defined, combinable instruments.

All instruments of the instrument set are intended for transient use.

## 6. Indications

### Implants

General Indications:

- Mobility-limiting diseases, fractures or defects of the hip joint or proximal femur which cannot be treated by conservative or osteosynthetic procedures

Indications:

- Primary and secondary osteoarthritis
- Rheumatoid arthritis
- Correction of functional deformities
- Avascular necrosis
- Femoral neck fractures
- Revision after implant loosening dependent on bone mass and quality

In addition to the indications and general indications for the Multi Hole Shell

- Where primary stability can only be achieved with the positioning of additional bone screws

In addition to the indications and general indications for the Dual Mobility

Inserts & Shell/Insert Adapters

- Dislocation risk

### Instruments

- The MobileLink Acetabular Cup System instrument set arrangement in the present combination is intended for the application and implantation of the MobileLink Acetabular Cup System.

## 7. Contraindications

### Implants

- Acute and chronic infections, local and systemic, insofar as they may compromise the successful implantation
- Allergies to (implant) materials
- Insufficient / inadequate bone mass- or quality, which prevents a stable anchorage of the prosthesis.

### Instruments

- The MobileLink Acetabular Cup System instrument set arrangement in the present combination is not intended for the application and implantation of another Medical Device aside from the MobileLink Acetabular Cup System.
- Material intolerance to the instrument material
- Instruments are not allowed to be implanted.

## 8. Possible Risks and Side Effects

- Fracture (bone, periprosthetic)
- Implant damage, implant fracture
- Infection
- Inflammation
- Instability, Dislocation
- Metal allergy
- Migration
- Malalignment
- Noises (e.g. cracking, popping, snapping, squeaking or grinding)
- Pelvic discontinuity
- Prosthetic impingement
- Septic, aseptic loosening
- Soft tissue problems
- Wear
- Residual complaints

## 9. Clinical Benefit

The clinical benefit of the MobileLink Acetabular Cup System is defined as

- Relieve of pain, compared to untreated diseased hip joint
- Increased joint functionality, compared to untreated diseased hip joint
- Increased stability of the joint, compared to untreated diseased hip joint
- Increased joint mobility, compared to untreated diseased hip joint

The clinical benefit of the instrument set arrangement for the MobileLink Acetabular Cup System is the implantation of the MobileLink Acetabular Cup System.

## 10. Implant Material

Please refer to the relevant surgical technique associated with the system and the identification on the packaging for further information on implant materials. Further information on the material compositions is available from the manufacturer upon request.

- Cobalt-based alloy, CoCrMo, ISO 5832-12; ASTM F1537, (EndoDur-S)
- Titanium-aluminum alloy, Ti6Al4V, DIN EN ISO 5832-3 ;ASTM F136, (Tilastan-S)
- Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene, UHMWPE, ISO 5834-2 / ASTM F-648 / ASTM F-2565, (X-LINKed PE)
- Highly cross-linked polyethylene based on ultra-high molecular weight polyethylene, UHMWPE, with vitamin E, ISO 5834-2 / ASTM F-648 / ASTM F-2565 / ASTM F-2695, (E-DUR)
- Calcium phosphate coating, ASTM F-1609, (HX)
- Titanium-aluminum alloy, Ti6Al4V, DIN EN ISO 5832-3, (Tilastan-E)
- Zirconia Toughened Aluminium Oxide, ZTA, ISO 6474-2
- Commercially Pure Titanium coating, CpTi, DIN ISO 5832-2, ASTM F1580

## 11. Use of Ceramic Inserts

- Ceramic inserts sterilized by gamma rays may have changed colour. This has no effect on the strength or any other property of the ceramic inserts.
- The ceramic insert is extremely susceptible to damage. Even small scratches or points of impact may cause excessive wear or breakage and give rise to complications. Extremely careful handling is therefore advised. Never strike the ceramic insert directly with a metal hammer.
- For acetabular shells in retroversion, a ceramic insert must not be used. Otherwise, increased surface pressure may develop on the edge of the cup and thus to grain break-out from the ceramic insert associated with increased ceramic debris. Excess ceramic debris can cause tissue reactions, prosthesis loosening, and, in extreme cases, to breakage of the ceramic insert.
- Ensure adequate joint tension is achieved on implantation.
- Before inserting the ceramic insert, the inside of the acetabular shell must be carefully cleaned and checked to ensure that the surrounding tissue does not impede insertion of the ceramic insert. To fix the ceramic insert in the acetabular shell, please note the information in the surgical technique for the combined implant system.
- On very rare occasions, in vivo fracturing of the ceramic insert may occur.
- If a ceramic component breaks, a pairing of metal (ball head) with polyethylene (insert) or of metal with metal must not be used for revision.

## 12. Instrument Material

Please refer to the REF list in this document and to the relevant surgical technique associated with the system and the identification on the packaging for further information on instrument materials.

Further information on the material compositions is available from the manufacturer upon request.

## 13. CMR Substances

Some system components contain cobalt as an alloy ingredient in a concentration above 0.1 % weight by weight.

Cobalt is listed as a CMR (carcinogenic, mutagenic and toxic to reproduction) substance.

The hazard class and category code(s) for cobalt are:

- Carc. 1B
- Repr. 1B

For identification of the affected components and further material specification, please refer to the REF list and the legend of materials in this document.

## 14. Implant Selection, permissible Combinations

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging.

Also, refer to the following sections of this document and to the relevant surgical technique associated with the system for further information on implant selection and permissible combination and for information on allocating and handling the instruments to be used for the implantation.

Combinations with implants from other manufacturers and / or combinations with LINK implants that deviate from the surgical technique specifications have not been tested and are not permitted.

Ceramic inserts may only be combined with the prosthesis components that are listed in the associated surgical technique of the combined implant system. The wear couple consists of two articulating joint surfaces of precisely defined geometry and precisely defined material. Combining BIOLOX components with LINK CeraDur components is not permitted.

## 15. Permissible materials for tribological pairings of implants

The implant components of the system only form tribological pairings when combined with other Link implants, e. g. prosthesis heads.

When combined with other Link implants, the following tribological pairings may be formed:

- CoCrMo alloy/X-LINKed or E-Dur
- Composite ceramics/X-LINKed or E-Dur
- Composite ceramics/ composite ceramics

“Ceramic on metal” and “metal on metal” pairings are prohibited within the area of the bearing surface. The same applies for the pairing of BIOLOX delta with a ceramic component by another manufacturer.

## 16. Range of Motion of the Implants

The range of motion of the standard combinations corresponds to DIN EN ISO 21535.

The range of motion rises as the head diameter of the prosthesis increases.

To reduce the risk of dislocation in patients with an increased tendency to luxation, Waldemar Link GmbH & Co. KG offers components that inhibit luxation (shouldered, anti-dislocation, inclined). These include, for example, Cup Inserts with shoulders (raising the cup edge) and Dual Mobility Inserts.

There is less range of motion than with a standard combination because the neck of the head and the shoulder reduce the range of motion.

The positioning of the prosthesis components has a direct influence on the range of movement and thus potentially on the risk of impingement, luxation, or subluxation. If the acetabular shell inclination is too steep, this increases the load on the edge of the cup. Increased wear and tear may develop as a result.

The inclination of the cup components should not be significantly above or below 40°–45°.

The anteversion of the cup components should not be significantly above or below 10°–20°.

The range of motion will be reduced outside these limits. This can lead to subluxations and/or dislocations of the ball head out of the ceramic insert.

## 17. Implant Diameters

Please refer to the relevant surgical technique associated with the system for further information on implant diameters.

## 18. Implant Anchoring

The MobileLink Acetabular Cup System is implanted without cement.

## 19. Lifetime

The lifetime of our implants is limited in principle and is determined by individual factors such as, for example, body weight and the level of activity of the patient, as well as by the quality and professional execution of the implantation. Based on these individual influencing factors, Waldemar Link defines the overall average lifetime of an implant based on its survival rate (i.e. the proportion of functional implants after a certain period of time starting from the time of implantation). According to the results of the tests performed, the survival rate of our implants corresponds to the general state of the art at the time of approval of the implants.

For the lifetime of the instruments please refer to the description in the reprocessing instructions H50.

## 20. Reprocessing / Reuse

The implants are supplied as sterile single-use devices. Implants whose protective packaging is opened or damaged, or implants which have already been implanted are not permitted to be reprocessed or reused. Unforeseeable changes may occur in ceramic implants during sterilization.

Likewise, a damaged ceramic insert must not be used but must be disposed of. This also applies, for example, to a ceramic insert that has fallen on the floor.

If ceramic inserts have already been used, there is a risk that they may have suffered damage that may not be visible to the naked eye. Damage of any kind can negatively affect the functionality and/or stability of the ceramic. Safe use can therefore not be guaranteed.

Therefore, only originally packaged, unused, and undamaged ceramic inserts may be implanted.

The following risks may occur if implants are re-used:

- Infections
- Reduced implant lifetime
- Increased wear and wear debris complications
- Disease transmission
- Inadequate implant fixation
- Limited implant function
- Implant response and / or rejection

Instruments must be disinfected and sterilized prior to use. For more information please refer to the related chapters in this document and to the description in the reprocessing instructions H50.

Additionally, observe our separate packaging and cleaning instructions for instruments.

Single-use products may not be reused.

## 21. Resterilization

The Implants are designed for single-use only. Resterilization is not permitted. Implants, as well as their materials are not suitable to be resterilized.

Unpredictable degradations may occur in these implants during resterilization. For sterilization of the Instruments, please refer to the description in the reprocessing instructions H50.

## 22. Storage and Transportation

Sterile-packaged implants must be stored in the undamaged original packaging in buildings with adequate protection against damage due to impacts, frost, humidity, excessive heat, and direct sunlight.

For Storage and Transport of the Instruments, please refer to the description in the reprocessing instructions H50.

## 23. Information for Patient Advisory

If the implantation of this system is considered to be the best solution for the patient and one of the circumstances described in section 24 is applicable to the patient, it is necessary to advise the patient with regard to the anticipated effects that these circumstances could have on the success of the operation. It is further recommended that the patient be informed about measures that he or she can take to reduce the effects of such complications. All information provided to the patient should be documented in writing by the operating surgeon. An implant ID for the patient must be handed over by the surgeon or hospital and the patient must be informed about the availability of a special patient information.

The patients should also be instructed:

- in detail about the surgery-related risks.
- in detail about the limitations of the implants, especially about the effects of excessive stress caused by body weight and physical activity, among other things. They should be encouraged to adjust their activities accordingly.
- about possible postoperative complications.
- about the material composition of the implant.
- that implants may respond to metal detectors during security checks (e.g. at airports) and carrying an implant ID as a proof is recommended.
- that implants may interact with medical imaging technique (e. g. MRI)

## 24. Circumstances that can interfere with the Success of an Operation

- Severe osteoporosis
- Severe deformities
- Local bone tumors
- Systemic diseases
- Metabolic disorders
- Case history of infections and falls
- Drug dependency or abuse, including excessive alcohol and nicotine consumption
- Obesity
- Mental disorders or neuromuscular diseases
- Heavy physical activities associated with strong vibrations
- Hypersensitivities

## 25. Warnings / Precautions

- The reuse of LINK single-use products is not permissible.
- Implants must be handled with great care and should not be modified or changed, even the smallest scratches and damages can considerably impair their stability or performance. Damaged implants are not permitted to be used.
- Surfaces provided for the connection of modular prosthetic components (cone, pins, screws) must not be damaged and may need to be cleaned with sterile liquid and dried before being joined together, so that neither blood nor any other coating impairs any of the connections, which could compromise the reliability of the connection.
- Do not manipulate or misuse instruments. We do not accept liability for products that have been modified, subjected to unintended use, or used improperly.
- For the processing of LINK instruments, it is presumed that the personnel have technical knowledge level I (Germany) and in other countries technical knowledge and expertise.
- Medical devices that are sent in for servicing must be processed beforehand in such a way that they cannot constitute a hazard to third parties.
- Products made of plastic (e.g. polyamide (PA), polyethylene (PE), polyoxymethylene (POM), ultra-high molecular weight polyethylene (UHMWPE)) may not be localizable using external imaging procedures.

## 26. Preoperative Planning

Preoperative planning provides important information to identify the appropriate implant system and select the components of a system. Make sure that all components required for the operation are laid out and ready in the operating room. Trial implants to verify proper fit (where applicable) and additional implants should be kept at the ready, in case other sizes are needed or the intended implant cannot be used. All LINK instruments necessary for the implantation must be on-hand, sterilized and intact.

If prosthesis implantation is indicated, then it must be taken into consideration, along with the overall circumstances of the patient:

- that all non-surgical and surgical treatment alternatives for the joint disease have been considered
- that artificial joint replacement performance is categorically inferior to natural joint performance, and an indication-related improvement in the preoperative condition is the only aim here
- that proper selection, placement and fixation of the devices are decisive factors, which will determine the life of the implant.
- that an artificial joint may loosen due to stress, wear and tear, and infection, or luxation or dislocation may occur
- that revision surgery, which under certain circumstances may exclude the possibility of restoring joint function, may be necessary due to loosening of the implant
- that the patient consents to undergo the operation and accepts the risks involved
- that if load-transferring bone structures are damaged, then the loosening of the components, bone and implant fractures, as well as other serious complications cannot be ruled out
- that if the patient is suspected of having allergies and tests positive on the applicable tests, then the patient's foreign body sensitivities (material tolerances) must be examined
- that acute and chronic infections, local and systemic, may compromise the successful implantation, so pre-operative microbiological analysis is recommended

Generally, the mechanical failure or fracture of an implant is a rare exception. However, this cannot be excluded with absolute certainty despite the sound structure of the implant.

This may be due to stress on the implant and prosthesis as the result of a fall or accident, among other things.

If the bone area where the implant is anchored is altered in such a way that the prosthesis is no longer able to withstand normal stress and an area of the prosthesis becomes subject to a stress imbalance, then a mechanical failure of the implant system may result. Such stress imbalances may also occur if the anchoring elements of the implants are required to form a bridge over larger bone deficiencies without optimal reinforcement of the bone. It is recommended that the implant with the largest possible anchoring elements be used. Proper preparation for surgical procedures also includes the functional testing of implants and instruments prior to use.

For definitive identification information on the product such as system compatibility, article number, material and shelf life, refer to the identification on the implant and / or the packaging. You should also take advantage of the training courses and printed materials provided for your information. To learn more, please contact the Waldemar Link GmbH & Co. KG sales office or your field representative.

## 27. Handling

All implant components are supplied sterile as single-use devices in individual packages. The implant components are sterilized by gamma sterilization, at least 25 kGy.

By contrast, components made of highly cross-linked polyethylene or highly cross-linked polyethylene with vitamin E are sterilized with ethylene oxide (ETO).

Implants should always be stored in their unopened protective packaging. Examine the packaging for damage before using the implant. Damaged packaging can have an adverse effect on both the sterility of the device as well as the proper performance of the implant, such that the device may no longer be used.

- Check the use by date on the implants. Implants with expired use by dates are no longer permitted to be used for implantation!
- After opening the package, check that model and size of the implant are matching with the information printed on the package labeling.
- Observe the pertinent standards for the aseptic handling of devices during and after removal of the implant from the packaging.
- When removing the packaging, make a record of the batch or serial numbers on the label, since this information is decisive for batch tracing. Self-adhesive labels with this information are enclosed with every package for your convenience.

Instruments must always be treated with care, this particularly applies during transport, cleaning, maintenance, sterilization, and storage. The sterile status of the instruments depends, inter alia, on the sterile items packaging and the prevailing storage conditions and must be established together with the operator's hygiene officer on a case-by-case basis. Direct sunlight must be avoided. Improper handling and care, as well as unintended use, can lead to premature wear or damage.

Devices made of plastics (e.g. PP-H) may not be located by means of an external imaging device.

## 28. Intraoperative Use

Please refer to the relevant surgical technique associated with the system for information about the intraoperative use of the system.

## 29. Postoperative

In addition to movement and muscle training, special attention must be paid to carefully instructing the patient during the postoperative phase.

Physician-supervised postoperative monitoring of healing progress is recommended. Where applicable, patients should also be advised on how to avoid overstraining themselves.

Follow-up examinations should be carried out regularly or immediately if symptoms occur.

## 30. Notes on MRI and CT Examination Procedures

Our implants were not evaluated for safety and compatibility for MRI and CT examination procedures.

In the case of our metallic implants and implant components, MRI examinations pose potential risks to the patient due to possible heating and migration of the implants or implant components.

Likewise, there is a potential risk of artefact formation in MRI and CT examinations of our metallic implants and implant components.

The probability of occurrence and the extent of the potential risks mentioned depend on the type of device used, its device parameters and the sequences used.

Always follow the instructions in the operating instructions of the manufacturer of the device used for the imaging.

The selection of the imaging examination procedure and the assessment of possible side effects is the responsibility of the examining physician.

The examining physician must consider the individual condition of the patient and other diagnostic methods.

## 31. Explantation of Implants / Revision Surgery

Please refer to the relevant surgical technique associated with the system for information about explantation of the implants and revision surgery.

## 32. Disposal

Packaging and system components to be discarded must be handled in compliance with your national and local regulations for hospital disposal.

## 33. Instruments

Please refer to the description in the reprocessing instructions H50 for the:

- initial use
- performance test
- maintenance
- manual cleaning
- cleaning in a washer-disinfector
- reprocessing
- sterilization
- servicing
- transport

## 34. Link to Summary of Safety and Clinical Performance

At the time of creation of this document, the EUDAMED database was not yet active. Therefore, no link to the summary of safety and clinical performance can be given here.

## 35. Requests

Requests of any kind should be directed to Waldemar Link GmbH & Co. KG (see contact information in this document).

## 36. Complaints about our Products

All complaints must be addressed to Waldemar Link GmbH & Co. KG at: [complaint@link-ortho.com](mailto:complaint@link-ortho.com).

In the event of a complaint, the name or reference number of the corresponding component should be specified with the serial number (SN) or the lot number (LOT), your name, and your contact address. The reason for the complaint should be given in brief.

## 37. Report of Serious Incidents

Any serious incident that occurs in relation to the device must be reported to the manufacturer and the authority responsible for your location.

Waldemar Link GmbH & Co. KG and/or other corporate affiliated entities own, use or have applied for the following trademarks in many jurisdictions: LINK, BiMobile, SP II, Modell Lubinus, E-Dur, EndoDur, T.O.P. II, BetaCup, CombiCup PF, CombiCup SC, CombiCup R, MobileLink, C.F.P., LCU, SP-CL, LCP, MIT-H, Endo-Model, Endo-Model SL, MP, MEGASYSTEM-C, GEMINI SL, SPAR-K, LCK, Link OptiStem, HX, TiCaP, X-LINKed, PorAg, PoroLink, LINK PorEx, BiPorEx, PorEx-Z, SatinLink, TrabecuLink, Tilastan, customLINK, RescueSleeve, Stactip, VACUCAST.

Other trademarks and trade names may be used in this document to refer to either the entities claiming the marks and/or names or their products and are the property of their respective owners.

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